

## PA Affinity Columns

### Product Information

---

**Catalog Number**

ACPA0002

---

**Size** 8.9×63mm

---

**Description**

PA Affinity Columns adopt porous silica as media, which are chemically conjugated with recombinant protein A via hydrophilic coating.

---

**Phase**

CD PA

---

**Application**

Use for the titer determination and small-scale purification of antibodies.

---

**Packing Material**

Porous silica

---

**Particle Size (μm )**

30 μm

---

**Pore Size (Å)**

1000 Å

---

**ID×Length (mm)**

8.9×63mm

---

**pH stability**

2-13

---

**Recommended operating pressure**

5,000 psi

---

**The maximum operating temperature**

---

## PA Affinity Columns

40

---

**Recommended flow rate range**

1500 cm/hour

---

**Dynamic Binding Capacity (mg/ml)**

67 mg/ml at flow rate 250 cm/h

---

**Immobilized ligand**

Recombinant protein A

---

**Shipping solution**

0.1 M sodium phosphate, pH8

---