

Capto L in ReadyToProcess columns

Product Information

Cat#No# Capto-005P

Product Overview

These ReadyToProcess chromatography columns are prepacked with Capto L chromatography resin and are excellent for the capture of antibodies and antibody fragments. ReadyToProcess columns are validated high-performance bioprocessing columns that are supplied prepacked and ready for use.

Characteristic

High binding capacity and selectivity for a broad range of antibody fragments reduce process time and amount of resin required.

Enables a platform approach for purification of antibody fragments.

Increases productivity, process economy, and process design flexibility due to the rigid agarose matrix.

Low ligand leakage.

Applications

For manufacturing of biopharmaceuticals for clinical phase I and II studies. Depending on the scale of operations, they can also be used for full-scale manufacturing, as well as for preclinical studies.

Maximum operating pressure

2 bar (0.2 MPa, 29 psi)

Recommended column height

200 mm (7.87 in)

Storage

2 to 8°C, 20% Ethanol

Pack size

2.5 L

Dimensions

126 × 200 mm

Capto L in ReadyToProcess columns

Column i.d.

126 mm
