

Zearalenone purification Kit (Magnetic beads Method)

Product Information

Alias

Zearalenone; Purification Kit; Magnetic beads

productkeywords

Zearalenone; Purification Kit; Magnetic beads; Zearalenone purification Kit (Magnetic Beads Method); KIT-026

Product Overview

Mycotoxin purification kit (Magnetic Beads Method) is a rapid and high-throughput purification kit for mycotoxin testing in grain and oil, food, feed, dairy, Chinese medicine and other samples. This kind of product uses magnetic beads with surface modification of antibodies to specifically adsorb mycotoxins in the sample extraction solution, and can complete the whole process of impurity washing and purification and elution of the target material with automatic purifier.

Storage

Sttorage at 2-8°C, do not freeze.

Product features

Efficient and fast: single sample processing time is less than 30min High throughput: Different models can handle 1-10 or 1-24 samples at the same time Automation: One-click start, no manual operation, convenient and fast, reduce human error.

Tel: 1-631-559-9269 1-516-512-3133

Fax:1-631-938-8127