

Pack-3 C18 Columns

Product Information

	normation.
Catalog Nur	ıber
18030026	
Size	2.1×250
Description	
	adopts full-coverage bonding silica filler, which has extremely high stability. These columns are with highly(100%) aqueous mobile phases.
Phase	
Pack-3 C18	
Application	
-	rations of both polar and non-polar compounds, such as pharmaceuticals, peptides, amino acids and organic acids.
Packing Mat	erial
Spherical, hig	h purity (<10 ppm metals)
Particle Size	(μm)
5	
Pore Size (Å)
120	
ID×Length (nm)
2.1×250	
pH stability	
3-7.5	
Recommend	led operating pressure
under 3,000	osi

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA



Pack-3 C18 Columns

The maximum operating temperature	
60	
Phase structure	
Monomeric and fully endcapped	
Pore volume	
1.0 ml/g	
Surface area	
300 m²/g	
% Carbon	
0.17	
Mobile phase compatibility	
Compatible with 100% aqueous mobile phase	

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA