

Pack-3 C18 Columns

Product Information

i i oddot ii	inormation
Catalog Num	ber
18030048	
Size	10.0×150
Description	
	dopts full-coverage bonding silica filler, which has extremely high stability. These columns are vith highly(100%) aqueous mobile phases.
Phase	
Pack-3 C18	
Application	
-	ations of both polar and non-polar compounds, such as pharmaceuticals, peptides, amino acident nd organic acids.
Packing Mat	erial
Spherical, hig	h purity (<10 ppm metals)
Particle Size	(μm)
5	
Pore Size (A	
200	
ID×Length (r	ım)
10.0×150	
pH stability	
3-7.5	
Recommend	ed operating pressure
under 3,000 p	si

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA



Pack-3 C18 Columns

The maximum operating temperature		
60		
Phase structure		
Monomeric and fully endcapped		
Pore volume		
1.0 ml/g		
Surface area		
200 m²/g		
% Carbon		
0.1		
Mobile phase compatibility		
Compatible with 100% aqueous mobile phase		

Fax:1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA