

## Pack-1 C8 Columns

## **Product Information**

Catalog Numb	ber
C8010009	
Size	4.6×150
Description	
Pack-1 C8 has wide hydrophol	lower hydrophobicity and higher efficiency. Suitable for the separation of compounds with a bic range.
Phase	
Pack-1 C8	
Application	
Use for separa organic acids.	tions of acidic, neutral and basic organic compounds, such as pharmaceuticals, peptides, and
Packing Mater	rial
Spherical, high	purity (<10 ppm metals)
Particle Size (	μm )
3	
Pore Size (Å)	
120	
ID×Length (m	m)
4.6×150	
pH stability	
3-7.5	
Recommende	d operating pressure
under 3,000 ps	i

Tel: 1-631-559-9269 1-516-512-3133

Email:info@creative-biomart.com

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA



## Pack-1 C8 Columns

The maximum operating temperature		
60		
Phase structure		
Monomeric and fully endcapped		
Pore volume		
1.0 ml/g		
Surface area		
300 m²/g		
% Carbon		
0.11		

Tel: 1-631-559-9269 1-516-512-3133

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA