

## Streptavidin HP SpinTrap Buffer Kit

### Product Information

<b>Cat#No#</b>	St-382C
----------------	---------

### Product Overview

Pre-made buffers optimized for the affinity-based enrichment of a target protein/immunoprecipitation using Streptavidin HP SpinTrap columns.

High quality buffers.

All buffer components listed and defined.

Convenient format - stock solutions can be diluted directly into the provided bottles for storage.

Eliminates time-consuming buffer preparation, which allows fast, reproducible protein enrichment.

### Storage

4 - 30°C

### Binding buffer

0.5 M Tris and 1.5 M NaCl, pH 7.5

### Elution buffer

1 M Glycine-HCl, pH 2.9

### Wash buffer

0.5 M Tris and 1.5 M NaCl, pH 7.5