

Blood RNA Extraction Kit

Product Information

Alias

RNA; Blood RNA; Extraction Kit

productkeywords

RNA; Blood RNA; Extraction Kit; Blood RNA Extraction Kit; KIT-019

Product Overview

The Blood RNA Extraction Kit uses a new magnetic bead technology to quickly and efficiently extract highquality RNA from anticoagulant treated fresh whole blood samples. This product uses superparamagnetic microspheres, in a special buffer system, can specifically adsorb nucleic acids, separate and purify to obtain high quality RNA.

Storage

Sttorage at 2-8°C, do not freeze.

Product features

Full automation: One step nucleic acid extraction, can be fully automated operation.

Simplicity: No red blood cell cracking step required, no heating operation.

Wide applicability: Suitable for different anticoagulant treatment of blood samples.

High safety: no irritating odor reagents, no phenol chloroform and other organic solvents.

Tel: 1-631-559-9269 1-516-512-3133

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA